

TOUCH SCREEN MOISTURE BALANCE MA-SERIES



FEATURES

- Four languages menu (Spanish, Germany, English, Chinese)
- HBM Load cell Base;
- 5 inch touch panel;
- 40-180°C test temperature
- Built in RS232 for PC/printer or other equipment;
- High Precision, fast measure, stable accuracy, simple operation
- Three drying mode
- Save 32 groups data

Contact: sales@sassoscales.co.ke

Website: www.sassoscales.co.ke

Display 5inch Touch Screen Moisture Balance

Capacity: 60g

Readability: 0.001g

Moisture Range: 0.00% to 100.00%

Pan Size 90mm

Heat Source: Halogen Lamp

Temp. Range 40°C to 180 °C

Moisture Range 0.00%-100.00%

Time Setting 0-99 minutes, interval 1 minute

Number of Storage 32 groups

External

Calibration

Connectivity Configurable RS 232 interface, optional USB