



## PRECISION SCALES

### JA SERIES- CLASS III

MODEL	CAPACITY	READAILITY (g)	PAN SIZE
JA-220	220g	0.1mg	DIAL
JA-303H	310g	0.1mg	DIAL
JA-1000	1000g	0.001g	DIAL
JA-5000	5000g	0.001g	DIAL



JA-SERIS are micro precision digital weighing scales that are simple to operate, clean and maintain. They have a large bright LCD display comfortable for your readings and gives very stable readings to avoid errors. JA series are doing well in the market for its ability to do internal calibration thus maintaining the accuracy of the machine. Also with them are serial port connectors for printer and computer connection for easier data management. SUITABLE FOR LABORATORIES USES

### CONTACT US

#### Physical Address

Scale And Services Solutions  
 P.O BOX 11333, Nairobi  
 North Airport Road, Embakasi

Email: [support@sassoscales.co.ke](mailto:support@sassoscales.co.ke)

**Need Help / Quick Contacts**

+2547 24 895 295

[sales@sassoscales.co.ke](mailto:sales@sassoscales.co.ke)

[support@sassoscales.co.ke](mailto:support@sassoscales.co.ke)